Warmly welcome

- **Manufacturer of Animal Components Free recombinant proteins with *Eschericoli* system.**

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- **Supplier of recombinant protein to Biopharmaceutical company.**

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- Established in 2008.
- In 2009, the first product, recombinant carboxypeptidase B was successful in market.
- In 2010, both rCPB and rTrypsin gained full success in Chinese market.
exported rCPB and rTrypsin to Europe, India and Asia, got success in international market.
2010年·重组羧肽酶B和重组胰蛋白酶在国内外市场获得成功。
- In 2012, move to a new office, Juke Biology Park.
- In 2012, Protein A, got success in market.
2012年·重组蛋白A·在市场获得成功。
- In 2012, registred in NSF, ISO9001-2008 and fullfilled in Jul. 2013.
- 2012年·质量体系认证·2013获得认证证书·2014-2015·两次复审。
- In 2013, , RCPB and RPT, gained Shanghai hige-tech program .
- 2013年·RCPB and RPT申报并获得上海市高转项目。
- In 2014, Protein A, gained Shanghai hige-tech program
2014·SPA申报并获得上海市高转项目。
- In 2015, applied for Shanghai hige-tech incorporation
2015年·申报上海市高新技术企业·扩大生产。

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

YAXINBIO rTrypsin and rCPB

■ USE OF RECOMBINANT CARBOXYPEPTIDASE B AND RECOMBINANT TRYPSIN IN INSULIN MANUFACTURING

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
Products List

Cat. No.	Product Name	Function	Application	Charaters	Package
RCPB01	Recombinant Carboxypeptidase B	catalyzes lysine, arginine and histidine from C-terminal of polypeptides.	r-Insulin production, antibody C-terminal identification, sequencing, etc	Special Activity (unit/mg) NLT 170	10mg,100mg,1gr, or bulk.
RPT0201	Recombinant Trypsine (porcine)	endoproteinase, hydrolyzes polypeptides at the carboxyl side of lysine and arginine, comply with USP 2014	r-Insulin production, biopharmaceutical process, cell culture	Special Activity (USPU/mg) NLT 3800	10mg,100mg,1gr, or bulk.
SRT0202	Recombinant Trypsine Sequence	endoproteinase, hydrolyzes polypeptides at the carboxyl side of lysine and arginine	peptide mapping, fingerprinting, and sequence analysis	Special Activity (USPU/mg) NLT 6000	20µg, 100µg, 1mg
RHT03	Recombinant Trypsine (human)	endoproteinase, hydrolyzes polypeptides at the carboxyl side of lysine and arginine, comply with USP 2014	biopharmaceutical process, cell culture, cell dissociation, human cell therapy, stem cell, etc.	Special Activity (USPU/mg) NLT 3800	10mg,100mg,1gr, or bulk.
RTS04	Recombinant Trypsine Solution	endoproteinase, hydrolyzes polypeptides at the carboxyl side of lysine and arginine, comply with USP 2014	cell culture, cell dissociation, biopharmaceutical process, etc.	Activity 2000 BAEE unit/ml	10ml,100ml,1L or bulk
RSPA05	Recombinant Protein A	binds to most human and mouse IgG subclasses. can be coupled to solid separation media to purify polyclonal or monoclonal IgG antibodies. can be coupled to a variety of molecules (such as fluorescent molecules, enzyme markers, biotin, colloidal gold and radioactive markers) used in antibody test in the process of Western-biot, ELISA or immunohistochemical tests, etc.		Purity NLT 95%	10mg,100mg,1gr, or bulk.
RLA06	Recombinant Lipase A	catalyzes the hydrolysis of fats and oils with excellent enantioselectivity.(1) Hydrolysis of trans-3-(4-methoxyphenyl) glycidic acid methyl ester [(±)-MPGM] (-)-MPGM.(2) Hydrolysis of (±)-naproxen methyl ester to produce (-)-naproxen.(3)Catalysis of ester substitution reaction, etc.		Activity 1000 unit/gr, 2000 unit/gr	1gr, 10gr, or bulk.
REK08	Enterokinase (EK)	cleaves lysine C-terminal preceded by four aspartic acids: Asp-Asp-Asp-Asp-Lys	delete extra N-terminal fusion protein to gain full recombinant protein	Activity (unit/µl) 10-50	100unit, 500unit, 1ku, or bulk.
RPK09	Recombinant proteinase K	endopeptidase. digest native proteins, thereby inactivating enzymes	used in the process of DNA extraction, etc.	Activity NLT 30	10mg,100mg,1gr, or bulk.
RCT10	Recombinant Chymotrypsine	endoproteinase, hydrolyzes polypeptides at the carboxyl side of aromatic amino acids: Tyr, Phe and Trp.	protein digests for peptide mapping or protein identification by peptide mass fingerprinting or MS/MS spectral matching.	Special Activity (unit/mg) NLT 1000	10mg,100mg or bulk.
V813	V8 (Endoproteinase Glu-C)	endoproteinase, cleaves peptide bonds C-terminal to Glu and Asp	Insulin analysis, peptide mapping, fingerprinting, and sequence analysis	Activity (unit/mg) 20	1 unit, 50 µg, 2mg.

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



A. Sketch of structure of Different kinds of antibodies

B. Specific binding between SPA and IgG Fc fraction


C. Solidated SPA bind IgG

Particle size comparison of Protein A affinity resins


Sepromax A40 Plus

PS-DVB ; 40um monodisperse


MabSelect SuRe

Agarose ; 40-130um Mean 85um


ProSep Ultra Plus

CPG ; 40-130um Mean 65um


POROS MabCapture A


PS-DVB ; 30-70um Mean 45um

YAXINBIO Principle of SPA affinity gel to purify antibody

YAXINBIO SPA affinity gel to cure immunity disease


In vitro, SPA affinity column absorb extra antibodies produced and existed in immunity patient blood, it is one kind of immunity therapy.

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
YAXINBIO R-Protein A production Process

Dynamic Binding Capacity

Pr-A affinity resin	Column size	Column volume
Sepromax® A40	ID 4.6×50 mm	0.83 ml
Sepromax® A40 Plus	ID 4.6×50 mm	0.83 ml
MabSelect SuRe®	HiTrap 7.0×25 mm	1.0 ml

Test conditions

System	AKTA® Purifier10
Binding buffer	20 mM PB , 0.15 M NaCl , pH 7.4
Elution buffer	0.1 M Gly-HCl , pH 2.5
CIP reagent	0.5 M NaOH
sample	1.0 mg/ml hlgG


Based on same matrix, the DBC curves of A40 and A40 Plus show nearly same variation trend. Contrary to MabSelect SuRe, the DBC of Sepromax A40 Plus doesn't decrease significantly with the increase of flow rate.

Dynamic Binding Capacity

Residence time vs. DBC

$$RT = \frac{H \times 60}{F}$$

Due to column size difference between Sepromax A40 Plus and MabSelect SuRe, we suggest to compare base on residence time using the formula $RT = H \times 60 / F$ where F is flow rate (ml/min), H is column height(cm), RT is residence time(min).



Example:
To obtain 40mg/ml DBC, RT of Sepromax A40 Plus is 1.2min, and MabSelect SuRe is 2.4min. Under the 10cm-column height, flow rate of Sepromax A40 Plus is 500cm/h, MabSelect SuRe is 250cm/h

Choose Sepromax A40 Plus, gain high efficiency

Recovery test

"Load-Eluent" Recovery

With polymer substrate, Sepromax A40 Plus showed excellent hydrophilic and very low non-specific interaction as polysaccharide media(MabSelect SuRe).


There's no significant difference on "load-elution" recovery rate between Sepromax A40 plus and MabSelect SuRe!

Alkaline stability study

Alkaline stability study

CIP condition: 0.5M NaOH, contact time 10 min each cycle, DBC test after every 20 cycles.

No difference!


Lifetime study

150-cycle lifetime test


A 150 loading and elution cycles study with CIP in every cycle revealed that Sepromax A40 Plus possessed excellent durability under harsh conditions.

步骤	溶液组成	体积或时间
平衡	20 mM PB , 0.15M NaCl , pH7.4	5CV
上样	1mg/ml hIgG+ BSA/Lysozyme	RT=0.6min , 50%DBC
冲洗	20 mM PB , 0.15M NaCl , pH7.4	5CV
洗脱	0.1M Gly-HCl , pH3.0	5CV
CIP	0.1M NaOH	4CV=15min

Overlay of chromatograms from runs in lifetime test


Relative DBC of Sepromax® A40 Plus during 150 lifetime cycles


After 150 runs , DBC showed modest 20% decrease

Overlay of chromatograms of DBC tests

Application


Purification of Human IgG from human serum

Purification of hIgG from Human Serum using Sepromax® A40 Plus.

An example of purification of human IgG is shown above. Human IgG was obtained from diluted human serum on a Sepromax® A40 column. Non-reduced and reduced SDS-PAGE gel electrophoresis results are shown right.

- Column: 4.6 ID×50 mm L
- Sample: Diluted Human Serum
- Loading buffer: 20 mM pH 7.5 PB with 0.15 M NaCl
- Elution: 0.1 M pH2.5 Glycine HCl

- Lane:
1. Marker
 2. Diluted human serum
 3. Flow-through
 4. Elution

HCP removal efficiency test

Impurity removal test – host cell protein (HCP)

➤ **Material and Method**

Host protein CHO testing kit, manufacturer: Cygnus, Cat. No.: CM015

ELISA method : put standards、 samples、 selected impurities、 internal standard into marked holes (96 holes) to get incubated. Then put anti-CHO protein antibodies to go through color development. Finally use ELISA reader determine 450 nm light absorption values.

➤ **Results**

Generate the standard curve, calculate the internal standard recovery rate, and then determine the HCP content using the following equation:

$$\text{HCP content}(\text{ng/mg}) = R_{TS} / (C_{ST} / d_{TS}) \quad \text{或} \quad \text{HCP}(\%) = [R_{TS} \times 10^{-6} / (C_{ST} / d_{TS})] \times 100$$

Where : R_{TS} is the detected ELISA sample HCP conc. (ng/mL) , C_{ST} sample conc. (mg/mL) , d_{TS} the dilution factor

Removal test of HCP



❖ Impurity removal: test- host cell proteins (HCP)

HCP (ng/mg)	Sepromax A40 Plus	MabSelect SuRe
supernatant	1754.3	1716.5
elution	1.9	4.5
Purification factor	9.2×10^2	3.8×10^2

Conclusion: Sepromax A40 Plus is a little bit better than MabSelect SuRe at the removal of residual host cell proteins.

DNA removal efficiency test

❖ Impurity removal: test- host DNA

rDNA (pg/mg)	Sepromax A40 Plus	MabSelect SuRe
supernatant	1.5 × 10 ⁷	1.6 × 10 ⁷
elution	154.4	431.3
DNA Purification factor	9.7 × 10 ⁴	3.7 × 10 ⁴

Figure on the right shows that both Sepromax A40 rPA material and product from company A possess ability to remove DNA. The two lines showed identical trend due to large lowering factors. Yet, Figure on the right shows that DNA lowering factor before and after the affinity chromatography. Sepromax A40 rPA material has higher ability in removing CHO cell DNA than the product from company A.

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NSF Int

789 North Dixboro Roac

Authorized Registration and
*Company is audited for conform

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Carl Blazik,
Vice President, Technical
Operations, NSF-ISR, Ltd.

NSF International Strategic Registrations

789 North Dixboro Road, Ann Arbor, Michigan 48105 | (888) NSF-9000 | www.nsf-isr.org

Authorized Registration and /or Accreditation Marks. This certificate is property of NSF-ISR and must be returned upon request.
*Company is audited for conformance at regular intervals. To verify registrations call (888) NSF-9000 or visit our web site at www.nsf-isr.org

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