Scientific questions of peptide mapping analysis of protein/antibody 蛋白及抗体肽图分析中的科学问题

Suxia Li

Shanghai Yaxin Biotechnology Ltd Co.

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d) Peptide map

Selective fragmentation of the product into discrete peptides is performed using suitable enzymes or chemicals and the resulting peptide fragments are analyzed by HPLC or other appropriate analytical procedure. The peptide fragments should be identified to the extent possible using techniques such as amino acid compositional analysis, N-terminal sequencing, or mass spectrometry. Peptide mapping of the drug substance or drug product using an appropriately validated procedure is a method that is frequently used to confirm desired product structure for lot release purposes.

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a) Peptide map

Selective fragmentation of the product into discrete peptides is performed by using suitable enzymes or chemicals. The resulting peptide fragments are analysed by high-performance liquid chromatography (HPLC) or other appropriate analytical procedures.

Page 77

The peptide fragments should be identified as far as possible using appropriate techniques such as mass spectrometry (MS) methods (e.g. electrospray ionization MS, matrix-assisted laser-desorption ionization time-of-flight MS). The use of MS/MS coupling should also be considered as it could reveal more detailed sequence information about the analysed peptide fragment. If one fragmentation method does not deliver the complete amino acid sequence, the use of an orthogonal enzyme or chemical cleavage method can increase the sequence coverage. The correct formation of the disulfide bridges may be characterized by the use of peptide mapping under reducing and non-reducing conditions.

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Journal of Pharmaceutical and Biomedical Analysis
21 (2000) 1099–1128

JOURNAL OF
PHARMACEUTICAL
AND BIOMEDICAL
ANALYSIS

www.elsevier.com/locate/jpba

Validation of a peptide mapping method for a therapeutic monoclonal antibody: what could we possibly learn about a method we have run 100 times?☆

Jacob Bongers ^{a,*}, John J. Cummings ^a, Mary Beth Ebert ^a,
M. Marcia Federici ^c, Linden Gledhill ^a, Dipti Gulati ^a, George M. Hilliard ^a,
Brian H. Jones ^a, Kwan R. Lee ^d, Jacek Mozdzanowski ^a, Michael Naimoli ^b,
Sudhir Burman ^a

^a Department of Analytical Sciences, SmithKline Beecham Pharmaceuticals, 709 Swedeland Road, King of Prussia,
PA 19406-0939, USA

^b Department of Biopharmaceutical Production, SmithKline Beecham Pharmaceuticals, 709 Swedeland Road, King of Prussia,
PA 19406-0939, USA

^c Department of Pharmaceutical Technologies, SmithKline Beecham Pharmaceuticals, 709 Swedeland Road, King of Prussia,
PA 19406-0939, USA

^d Department of Statistical Sciences, SmithKline Beecham Pharmaceuticals, 709 Swedeland Road, King of Prussia,
PA 19406-0939, USA

Accepted 2 August 1999

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Journal of Pharmaceutical and Biomedical Analysis
16 (1997) 631–640

JOURNAL OF
PHARMACEUTICAL
AND BIOMEDICAL
ANALYSIS

Rapid analytical monoclonal an

K. Kannan ^a, M.C.
K. Seweri

^a Department of Quality Control C

^b Department of Analytical Ch

Rec

3.1. Protein unfolding, reduction, S-carboxymethylation and trypsin digestion

A concentration of 6 M guanidine hydrochloride was used for unfolding the monoclonal antibody. The extent of unfolding was monitored by circular dichroism and was found to be complete under the experimental conditions. The use of another common protein unfolding reagent, 8 M urea, resulted in incomplete digestion, most likely due to incomplete unfolding of the antibody.

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BIO PROCESS TECHNICAL

Solution-Phase Sample Preparation Approach for Peptide Mapping of Biological Therapeutics

Vajira Nanayakkara, William E. Werner, and Thomas Theriault

Table 2: Summary of data for protein 1, comparing the PD10 column buffer exchange, dialysis, and ultrafiltration methods

	PD10	Dialysis	Ultrafiltration
Sequence coverage by AA residues (%)	75	80	94
H chain C — terminal Lys clipped	Yes	Yes	Yes
H chain Asn 34, 99, and 349 glycosylation	Weak signal	Weak signal	Yes
H chain N — terminal deamidation	Not found	Weak signal	Yes
Salt adduct ions	Yes	Yes	None or trace
Data quality and reproducibility	Poor	Poor	Excellent
Sample processing time per sample	36 hours	36 hours	12 hours

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Rituximab USP Medicine Compendium Final Authorized Version 1.0 (2012)

Trastuzumab USP Medicine Compendium Final Authorized Version 1.0 (2013)

Bevacizumab USP Medicine Compendium Final Authorized Version 1.0 (2014)

Published on USP Medicines Compendium (<https://mc.usp.org>)

Rituximab

Final Authorized Version 1.0

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• B. PEPTIDE MAPPING

Use a chromatographic system. (Proceed as directed in *Biotechnology Derived Articles—Peptide Mapping <1055>*.)

Analyze the material to be tested by a chromatographic technique capable of resolving peptides generated from a Trypsin digest.
The digest is carried out under reducing conditions which provides NLT 80% digestion. The test procedure used provides a minimum of 90% coverage of the protein sequence.

Standard solution: Digest and dilute a portion of USP Rituximab RS in an appropriate diluent.

Sample solution: Digest and dilute a quantity of Rituximab in an appropriate diluent to obtain a nominal concentration of Rituximab similar to that of the *Standard solution*.

Control solution: Digest and dilute a portion of an appropriate control (non-Rituximab monoclonal antibody) in an appropriate diluent to obtain a nominal concentration of the control that is similar to that of *Standard solution*. [Note—The digests described in the *Standard solution*, *Sample solution*, and *Control solution* are conducted at the same time, using the same stock and concentration of reagents.]

IDENTIFICATION

• A. PEPTIDE MAPPING

Solution A: 0.1% Trifluoroacetic acid in water

Solution B: 0.1% Trifluoroacetic acid in acetonitrile

Solution C: 0.5 M dithiothreitol in water

Solution D: 0.5 M iodoacetamide in water

Solution E: 0.25 M tris buffer in water. Adjust with dilute hydrochloric acid to a pH of 7.5.

Solution F: 6 M guanidine hydrochloride and 1 mM EDTA in Solution E (denaturing buffer)

Solution G: 0.1 M tris buffer in water. Adjust with dilute hydrochloric acid to a pH of 7.8.

Solution H: 2 M urea in Solution G (digest buffer)

Solution I: 0.05 M acetic acid in water

Solution J: 1 mg/mL of trypsin in *Solution I*

Solution K: 10 mg/mL of USP Rituximab RS in water

Standard stock solution 1: Mix 100 µL of *Solution K*, 400 µL of *Solution F*, and 10 µL of *Solution C*, and incubate at 37° for 30 min. Add 24 µL of *Solution D* and incubate at room temperature for an additional 30 min in dark. Add 10 µL of *Solution C* and mix well.

Standard stock solution 2: Wash the PD-10 Sephadex G-25 column with 20 mL of water and equilibrate with 35–40 mL of *Solution H*. Load *Standard stock solution 1* on the column, and elute using *Solution H* in volumes of 700 µL each. Collect 6 independent fractions. Measure the absorbance of each fraction at 280 nm against *Solution H*. The fraction having an absorbance between 1.3 and 2.0 is used for digestion. [NOTE—If the absorbance of the fraction is more than 2.0, dilute it using *Solution H* to get an absorbance of 2.0.]

Standard solution: Mix well 50 µL of *Standard stock solution 2* and 2 µL of *Solution J* and incubate for 18–20 h at room temperature. Add 1 µL of trifluoroacetic acid and store the mixture at 4°.

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Questions:


- 1) Different buffer for denaturing (6M GuHCl) and for digest (2M urea)
- 2) Mode for buffer changing: G25 column
- 3) Digestion time: 18-20h in room temperature

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Rapigest SF surfactant from Waters

Similar as SDS, but more effective, and good activity for trypsin in Rapigest

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
Rapigest SF surfactant from Waters

Suggested procedure for in-solution digestions

1. Suspend lyophilized Rapigest SF powder in 1mL of 50mM NH₄HCO₃ to give 0.1%(w/v).
2. Suspend protein pellet in
3. Add DTT to 5mM
4. Heat the sample at 60°C for 30 min
5. Cool the sample
6. Add IAA to 15mM and place the sample in dark 30min
7. Add enzyme for digestion (1:100 to 1:20, w/w)
8. Incubate at 37 fro 1hr to overnight depending upon protein hydrophobicity

Sequencing grade trypsin:

TPCK-trypsin

Promega

Sigma

Recombinant trypsin for peptide map

YaxinBio

Roche

PRODUCT NAME	comments	Content of chymotrypsin in Trypsin
USP-trypsin	USP standard	<5%
Trypsin (sample)	Trypsin sample before TPCK treated	1.8%
TPCK-trypsin	Trypsin after TPCK treated	0.56%

There is chymotrypsin activity in TPCK-treated trypsin.

No any other enzymes activity.

So, it is unnecessary for TPCK treated.

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Page □18

Packaging

39(5) In-Process Revision: <89> ENZYMES USED AS ANCILLARY MATERIALS... Page 1 of 7

BRIEFING

< 89> Enzymes Used As Ancillary Material in Pharmaceutical Manufacturing. This new general test chapter provides analytical procedures to aid in the assessment of quality for enzymes that are used in biopharmaceutical manufacturing. Examples include trypsin, collagenase, pepsin, and papain. This chapter does not discuss the applications of these enzymes but rather focuses on tests to assess the qualities as process materials. Further, the chapter does not provide ways to limit the enzymes in the final medicinal product. The first enzyme discussed in the chapter is recombinant trypsin. Other enzymes will be added in subsequent revisions of the chapter.

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YAXINBIO Recombinant trypsin 重组胰蛋白酶

The activity of recombinant trypsin is determined using a chromogenic peptide substrate carbobenzoxy-valyl-glycyl-arginine-4-nitril-anilide acetate. The liquid chromatographic procedure in the test for *Purity* is based on analyses performed with the YMC-Pack ODS-A brand of L1 column. The retention time for the main peak is 12–17 min. A new Reference Standard, USP rTrypsin RS, is proposed for use during assessment of the system suitability for the proposed Assay and test for *Purity*.

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Comparation of trypsin

Trypsin


Test/Method	Specification 1:2500 in USP/CP	Specification 1:250
Appearance	White powder	Yellowish or brown
Solubility	Soluble	Soluble overnight
Microbial limits	-	none
<i>Staphylococcus aureus</i>	Negative	none
<i>Pseudomonas aeruginosa</i>	Negative	none
<i>Salmonella</i> species	Negative	none
Loss on drying	NMT 5.0%	none
Residue on ignition	NMT 2.5%	none
Limit of Chymotrypsin	NMT 5.0%	none
Activity	NLT 2,500 USP units/mg	NLT 250 USP units/mg

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
The retention time for the main peak for r-trypsin is 12-17min.
主峰的保留时间为： 12-17min

Retention time: NLT 1.0min
between two peaks of α -trypsin
and β -trypsin.

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NLT 70% for the peak area
of β -trypsin and NMT 20%
for the peak area of α -
trypsin.

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Activity of different concentration of rTrypsin in 50mM NH₄HCO₃, pH7.8, 37°C

ug/ml	10		20		50		100		200	
	U	100%	U	100%	U	100%	U	100%	U	100%
0h	110	100%	215	100%	600	100%	1140	100%	2300	100%
1h	115	105%	220	102%	575	96%	1170	103 %	2300	100%
2h	105	95%	240	112 %	540	90 %	1020	89%	2100	91 %
4h	120	109%	220	102 %	590	98 %	1160	102 %	1900	83 %
20h	75	68 %	135	63 %	290	48 %	440	39%	600	26%


Suggested used concentration: 10-50ug/ml

For example: 1ug in 50ul

Suggested digestion time: 4h

Stability of rTrypsin

**SDS-PAGE after incubated in 50mM NH₄HCO₃ pH7.8
37°C, 10h and 20h**

**No self-degradation. 无自降解
37°C, 50mM NH₄HCO₃ pH7.8 10h**

- 1: Marker
- 2: 10ug/ml 10h 0.2ug
- 3: 20ug/ml 10h 0.4ug
- 4: 50ug/ml 10h 1ug
- 5: 100ug/ml 10h 2ug
- 6: 200ug/ml 10h 4ug
- 7: 500ug/ml 10h 5ug
- 8: 10ug/ml 20h 0.2ug
- 9: 20ug/ml 20h 0.4ug
- 10: 50ug/ml 20h 1ug
- 11: 100ug/ml 20h 2ug
- 12: 200ug/ml 20h 4ug
- 13: 500ug/ml 20h 5ug

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YAXINBIO Stability of rTrypsin in GuHCl

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A组: rTrypsin activity (%) in 50mM NH₄HCO₃ +
2M GuHCl pH8.0, 37°C

ug/ml	10	20	50	100	200
0h	100	100	100	100	100
2h	100	100	98	95	100
4h	100	100	91	1000	100
6h	100	100	89	96	86
22h	100	93	68	90	71

D组: rTrypsin activity (%) in 50mM
NH₄HCO₃ + 1.5M GuHCl pH8.0, 37°C

ug/ml	10	20	50	100	200
0h	100	100	100	100	100
2h	100	100	100	100	100
4h	100	100	100	100	100
6h	100	100	100	100	95
22h	100	100	96	100	95

Stable in 1.5M or 2M GuHCl, pH8.0, 37°C

No self-degradation. 无自降解，活性稳定

Suggested process for pretreatment and trypsin digestion:

6M GuHCl in 50mM NH₄HCO₃ pH8.0 for denaturing and then dilute to 2M for trypsin digestion directly.

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





Yixin 5h

Promega 5h

■ Almost the same chromatography profile (TIC)

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
Enzyme: recombinant Carboxypeptidase B

ITEMS	SPECIFICATION
Source	Recombinant <i>E.Coli</i>
Purity	> 95%
Specific activity	>200 U/mg pro.
Trypsin content	<1ppm
Other enzymes	None
Form	Lyophilized powder
pI	5.4

STORAGE AND STABILITY

Recommend recombinant carboxypeptidase B lyophilized powder stored under 2-8°C in sealed container. It is stable within 24 months. After dissolved, it should be stored under -20°C, It is stable within 24 months and no activity lose after 10 times repeated freezing and thawing.

1, Merck CPB

2, YaxinBio rCPB

Enzyme: recombinant Carboxypeptidase B 重组羧肽酶B

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




Ratio of AB to CPB: 20:1 to 100:1

Used concentration: 0.5ug-1ug in 50ul

CPB peak due to
much more CPB
used (2.6ug in
loading sample).

YAXINBIO 测序级糜蛋白酶 r-chymotrypsin

Source	<i>E. Coli</i>
Purity	NLT 95% by HPLC analysis.
Physical form	Lyophilized powder
Specific activity	NLT 1500 USP units/mg pro.
No activity Contaminant	No carboxypeptidase A, or other proteases contaminant.

APPLICATION

Protein digests for peptide mapping applications or protein identification by peptide mass fingerprinting or MS/MS spectral matching. It is suitable for digestion reactions in-solution or in-gel.

RECOMMEND USAGE

To prepare 1-10mg/ml with 1mM HCl, used within 2 days, or stored below -20°C after repacked

STORAGE AND STABILITY

Recommend recombinant Chymotrypsin lyophilized powder stored under 2-8°C in sealed container. It is stable within 24 months. After dissolved, it should be stored under -20°C.

■ Sequencing grade enzymes for protein analysis

Cat. No.	Product Name	Function	Charaters	Package
SRCPB0102	Recombinant Carboxypeptidase B (Sequencing grade)	---↓Lys/Arg/His	Specific Activity (unit/mg)	100µg
			NLT 200	1mg
SRT0202	Recombinant Trypsin (Sequencing grade)	---Lys/Arg↓---	Specific Activity (BAEE unit/mg)	100µg
			NLT 15000 BAEE unit	1mg
SRCT10	Recombinant Chymotrypsin (Sequencing grade)	---Tyr/Phe/Trp↓--- at a lower rate, at Leu and Met.	Specific Activity (unit/mg)	100µg
			NLT 1500	1mg
V813	V8 (Endoproteinase Glu-C) (Sequencing grade)	---Glu/Asp↓---	Activity (unit/mg)	50µg (1U)
			20	2mg

YAXINBIO Animal free questions in vaccine

The first is : Non-serum media

The second is : recombinant trypsin


rTrypsin is to meet the needs of the cell culture and vaccine industry. With same performance to animal derived trypsin, rTrypsin provides efficient dissociation of cells from surfaces and tissues while maintaining cell viability and integrity.

Additional question is : recombinant trypsin inhibitor

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

Vaccination, early training, stick a
wall stage - **harvest cell** - culture

Trypsin is used during harvesting,
If propagate after harvesting, aprotinin
is needed to inhibit trypsin activity.

- ◆ Trypsin (EC3.4.21.4) has been available for many years as native enzyme isolated from porcine and bovine pancreas.
- 1:250 trypsin specially for cell dissociation
- 1:2500 trypsin collected by CP, USP, EP and JP etc in 1970s.
- 1:3800 recombinant trypsin issued by USP in 2014
- In 2015, CP2015 is preparing the standard for recombinant Trypsin.


HeLa Cell

YixinBio R-Trypsin

A1

A2


A3

B1

B2

B3

C1

C2

C3

A: rTryp1000 BAEE units/ml+0.01% EDTAA1: before; A2: after **3min** with rT;

A3: after 15h culturing

B: rTryp10000 BAEE units/ml+0.01% EDTAB1: before; B2: **after 30sec** with rT;

B3: after 15h culturing

C: rTryp2000 BAEE units/ml+0.01% EDTAC1: before; C2: after **2min** with rT;

C3: after 24h culturing

The effects of Recombinant trypsin components for cell detachment and its growth

Components	Time (min) for detachment	48 h Cell Growth condition
1000 BAEE units in PBS	7	+++
2000 BAEE units in PBS	6	+++
5000 BAEE units in PBS	4.5	+
10000 BAEE units in PBS	3	-
1000 BAEE units in PBS + 0.01% EDTA	3	+++++
2000 BAEE units in PBS + 0.01% EDTA	2	+++++
5000 BAEE units in PBS + 0.01% EDTA	1.8	+
10000 BAEE units in PBS + 0.01% EDTA	1.5	-

Recombinant Trypsin: YaxinBio
Cell Line: HeLa

1. The concentration of trypsin

If shift from 1:250 trypsin, convert with activity, such as,

if 0.25% for 1:250 (250 USP unit/mg),

to 0.025 for 1:2500 (2500 USP unit/mg),

to 0.016 for 1:3800 (3800 USP unit/mg).

2. The time for trypsin


If change FBS to non-serum medium, pay attention to the time, **to stop the trypsin with inhibitor.**

3. Different cell line owns different optimum concentration of trypsin, find a optimum low effective concentration


➤ Aprotinin /BPTI

While aprotinin and bovine pancreatic trypsin inhibitor (BPTI) are the same protein sequence, the term aprotinin is typically used when describing the protein derived from bovine lung.

58 amino acid, 3 disulfides

➤ recombinant Aprotinin /BPTI
--YaxinBio

Source	<i>E. Coli</i>
Purified by	HPLC
Physical form	Liquid in 0.1 M NaCl or White lyophilized powder
Specific activity	≥ 5 EPU/mg pro.
Purity	≥ 98% by SDS-PAGE
Contaminant activity	No any other protease contaminant.

Aprotinin-trypsin complex:
mol:mol

UNIT DEFINITION

One trypsin inhibitor unit (EPU) will decrease the activity of 2 trypsin units by 50% where one trypsin unit will hydrolyze 1.0 µmole of N-benzoyl-L-arginineethyl ether (BAEE) per sec at pH 7.6 at 25 °C.


A conversion factor for Aprotinin is: 1 EPU = 1 USP Aprotinin Unit = 1800 KIU.

Usage: equal- mol trypsin
(or 1/3 weight of trypsin)
store condition: 2-8°C or -20 °C

- soybean trypsin inhibitor
---from TrypZean sigma

Soybean trypsin inhibitor and other inhibitors work the same with TrypZean as they do with native trypsins (on a weight-to-weight basis)

Activity: >7,000 BAEE units/mg

Usage suggestion: 0.25-0.5 mg/ml in PBS.

Inhibit equal- volume trypsin

store condition: -5 °C -- -20 °C

180 amino acid, 2 disulfides

- **Manufacturer of Animal Components Free recombinant proteins with *Eschericholi* system.**

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- **Supplier of recombinant protein to Biopharmaceutical company.**

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Products List

Cat. No.	Product Name	Function	Application	Charaters	Package
RCPB01	Recombinant Carboxypeptidase B	catalyzes lysine, arginine and histidine from C-terminal of polypeptides.	r-Insulin production, antibody C-terminal identification, sequencing, etc	Special Activity (unit/mg) NLT 170	10mg,100mg,1gr, or bulk.
RPT0201	Recombinant Trypsine (porcine)	endoproteinase, hydrolyzes polypeptides at the carboxyl side of lysine and arginine, comply with USP 2014	r-Insulin production, biopharmaceutical process, cell culture		
SRT0202	Recombinant Trypsine Sequence	endoproteinase, hydrolyzes polypeptides at the carboxyl side of lysine and arginine	peptide mapping, fingerprinting, and sequence analysis	Special Activity (USPU/mg) NLT 6000	20μg, 100μg, 1mg
RHT03	Recombinant Trypsine (human)	endoproteinase, hydrolyzes polypeptides at the carboxyl side of lysine and arginine, comply with USP 2014	biopharmaceutical process, cell culture, cell dissociation, human cell therapy, stem cell, etc.		
RTS04	Recombinant Trypsine Solution	endoproteinase, hydrolyzes polypeptides at the carboxyl side of lysine and arginine, comply with USP 2014	cell culture, cell dissociation, biopharmaceutical process, etc.	Activity 2000 BAEE unit/ml	10ml,100ml,1L or bulk
RSPA05	Recombinant Protein A	binds to most human and mouse IgG subclasses. can be coupled to solid separation media to purify polyclonal or monoclonal IgG antibodies, can be coupled to a variety of molecules (such as fluorescent molecules, enzyme markers, biotin, colloidal gold and radioactive markers) used in antibody test in the process of Western-blot, ELISA or Immunohistochemical tests, etc.			
RLA06	Recombinant Lipase A	catalyzes the hydrolysis of fats and oils with excellent enantioselectivity.(1) Hydrolysis of trans-3-(4-methoxyphenyl) glycidic acid methyl ester [(±)-MPGM] (-)-MPGM.(2) Hydrolysis of (±)-naproxen methyl ester to produce (-)-naproxen.(3)Catalysis of ester substitution reaction, etc.		Purity NLT 95%	10mg,100mg,1gr, or bulk.
REK08	Enterokinase (EK)	cleaves lysine C-terminal preceded by four aspartic acids: Asp-Asp-Asp-Asp-Lys	delete extra N-terminal fusion protein to gain full recombinant protein		
RPK09	Recombinant proteinase K	endopeptidase, digest native proteins, thereby inactivating enzymes	used in the process of DNA extraction,etc.	Activity NLT 30	10mg,100mg,1gr, or bulk.
RCT10	Recombinant Chymotrypsine	endoproteinase, hydrolyzes polypeptides at the carboxyl side of aromatic amino acids: Tyr, Phe and Trp.	protein digests for peptide mapping or protein identification by peptide mass fingerprinting or MS/MS spectral matching.		
V813	V8 (Endoproteinase Glu-C)	endoproteinase, cleaves peptide bonds C-terminal to Glu and Asp	Insulin analysis, peptide mapping, fingerprinting, and sequence analysis	Activity (unit/mg) 20	1 unit, 50 μg, 2mg.

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- Manufacturer of AOF recombinant enzymes for recombinant human insulin

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